

Surface Adsorption of Chlorine Oxyanion Species to the Air-water Interface

Shane W. Devlin^{1,2}, Monika Blum^{2,3}


1: Department of Physics University of Nevada, Las Vegas, Las Vegas NV

2: Advanced Light Source, Lawrence Berkeley National Laboratory, Berkeley CA

3: Chemical Sciences Division, Lawrence Berkeley National Laboratory, Berkeley CA

Motivation

The adsorption of select ions to the liquid water interface has profound effects on important environmental processes, such as evaporation rates, aerosol chemistry, accelerated reactions, and more. Despite the overwhelming importance of these molecular details, the underlying physics of ion solvation models used to interpret surface phenomena are still widely debated. Chlorine oxyanions serve as a powerful model system, due to their abundance in seawater and importance in aerosol cloud chemistry. Here, to better understand how size, charge density, and hydration enthalpy, among other factors, influence surface adsorption, we compare the bulk and surface speciation of chlorine oxyanion species (ClO₄- and ClO₃-) in liquid water. Using Ambient Pressure X-ray Photoelectron Spectroscopy (APXPS) on planar liquid jets, we are able to track individual ionic species at various depths away from the interface.




All relevant species identified: O, Cl, and Na

Highly sensitive to concentration, even with

Ambient Pressure X-ray Photoelectron Spectroscopy of Liquid Flat Jets Colliding flat jet Hemispherical Analyzer Interaction Chamber Quartz nozzles Aperture to analyzer Liquid catcher Eiected photoelectron Temperature control of liquids Recirculating solution system Inner + outer capillaries for fast mixing and

Changing incident photon energy probes chemical species at different depths away from the interface


Characterizing Planar Liquid Jets Thin film interference models are used to reliably measure the flat jet thickness



kinetic studies

Better target for broad X-ray sources,

compared to cylindrical jets



200 eV

400 eV

600 eV

Larger KE probes deeper into the solution

0.5 M NaClO₄ and NaClO₃ solutions show similar surface speciation

300 400 500 600

On-going collaborative work: Soft X-ray Second Harmonic Generation at X-ray Free Electron Lasers

Conclusions

Combining surface selective measurements at XFELs with depth profiling at ALS for a "complete" description of ion distributions at the interface

Top of jet: ≈ 5 μm

Bottom of jet: ≈ 800 nm

0.5 M NaClO₄

0.5 M NaClO3



short acquisition times

